

Touchless Hand Hygiene Station

Dermados Station

Dermados Station Touchfree (5/20L)

Stand-alone unit with touch free operation for safe application of hand disinfectants.

.....

- For installation in entrance areas of office buildings, hotels, hospitals, food and beverage plants, life science and industrial customers, schools and kindergartens, supermarkets etc.
- Attracts people's attention for hand disinfection by hygienic High-End Stainless-Steel Design and clear guidance by pictograms.
- Maximum hand hygiene by touch free operation
- Battery driven system, no power cable required
- More than 6000 dosages possible (20 Liter jerrycan)
- Saves time and money as it is maintenance friendly with reduced refill checks


A user-friendly and safe Hygiene Station that can be easily installed wherever safe hand hygiene is required.



Dermados Station

Dermados Station Touchfree (5/20L)

- Convenient height for use
- Maximum hand hygiene by touch free operation
- Safe stand-alone unit due to battery driven system and no cable
- Economical use of 5/20 Liter Ecolab jerrycans by improved canister emptying
- Ensures high hygiene level by sufficient product supply of 5 or 20 L jerrycans
- Convenient option to lock single/multiple stations with existing key system
- Adjustable dosing rate (according hand disinfectant needs i.e. Hygienic Hand Disinfection acc. EN 1500)
- For frequent use areas
- Freely positional - no need for wall mounting or cable routing
 - ▶ Place it where it is needed!
- Easy maintenance
 - ▶ Easy standard battery exchange
 - ▶ Easy jerrycan exchange
 - ▶ Less inspections due to high product storage




TECHNICAL SPECIFICATIONS

Dermados Station	
EBS order no.:	10240138
Material	<ul style="list-style-type: none"> • brushed stainless steel • ABS plastic (dispenser) • silicone (tube)
Dosing rate	<ul style="list-style-type: none"> • 3 ml factory setting with LED double sensor technology • adjustable • one spare pump included
Battery power	4 x 1.5 V AA Mignon LR6
Battery duration	minimum 30.000 cycles with LED prewarning
Dimension (W x D x H)	370 x 617 x 1312 mm
Weight	approx. 20 kg

WE RECOMMEND THE USE OF

Skinman Soft Protect FF


Virucidal liquid hand disinfectant with Vitamin E, Glycerin & Panthenol



Hygienic Hand Disinfection acc. EN 1500 3 ml/20 sec.; EN 14476 (virucidal, incl. efficacy against Coronavirus and others) 30 sec.

Manodes GP

Liquid Hand Disinfectant for all areas except patient treatment and care in healthcare settings



Hygienic Hand Disinfection acc. EN 1500 2 x 3 ml/60 sec.; Effective against coronavirus 1 x 3 ml / 30 sec.

Use biocides safely. Always read the label and product information before use.

EMEA Headquarters
Ecolab Europe GmbH
 Richtistrasse 7
 8304 Wallisellen
 Switzerland
 Tel. +41 44 877 2000

Ecolab Engineering GmbH
 Raiffeisenstraße 7
 D-83313 Siegsdorf
 Tel. 02173 599 1888
 ee-service-de@ecolab.com

www.ecolab.com

ECOLAB®
 Everywhere It Matters.™