



AseptoStar



Powerful chlorinated alkaline cleaner and disinfectant for milking equipment and milk tanks.

- ▲ Powerful disinfection
- ▲ Strong performance in cleaning
- ▲ Best added value solution

APPLICATION

Pre-clean with water, clean and disinfect with **AseptoStar** at 0.4 - 0.7 % (minimum 10 min. at 40 - 50°C). Rinse with water of drinking quality. To be used in daily rotation with **HorolithStar**.



FEATURES AND BENEFITS

AseptoStar is a high-end alkaline chlorinated product with excellent cleaning and disinfection properties. It is suitable for all kinds of standard milking equipment, including cooling tanks. **AseptoStar** is in compliance with the Biocide Product Regulation (No. 528/2012).

The **Star** products are the top level program to ensure the best hygiene results.

SPECIFICATIONS

AseptoStar	
Physical form	Light yellow liquid
Active substance	Sodium hypochlorite
pH at 20°C:	13.5 - 14.0
Density at 20°C (g/cm ³):	1.18 - 1.22
Storage temperature:	-15°C to +30°C

LOGISTICS

Packaging	12 kg	24 kg	72 kg	245 kg
Product code	2321330	2321350	2334650	2321370
EAN code	4028162321333	4028162321357	4028162334654	4028162321371

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtstrasse 7
CH-8304 Wallisellen
www.ecolab.com