



# HorolithSolo

Powerful all-in-one cleaner for milking equipment and milk tanks.

- ▲ Powerful single-phase product
- ▲ Double concentrated powder
- ▲ Readily soluble

## APPLICATION

Pre-clean with water. Clean with **HorolithSolo** at 0.25% (minimum 10 min. at 40°C - 50°C). Rinse with water of drinking quality. Can be used for each cleaning.



## FEATURES AND BENEFITS

**HorolithSolo** is a strong stand-alone-cleaner which is suitable for all kind of milking equipment, including cooling tanks. It is able to eliminate all kinds of milk deposits. The double concentrated and readily soluble powder leads to less packaging waste and transport.

## SPECIFICATIONS

<b>HorolithSolo</b>	
Physical Form	White powder
Active Substance	Inorganic acid and blend of additives
pH at 20°C:	1.1
Density at 20°C (g/cm <sup>3</sup> ):	1.2
Storage temperature:	-10°C to +40°C

## LOGISTICS

Packaging	5 kg	15 kg
Product code	2321880	2321890
EAN code	4028162321883	4028162321890

A Material Safety Data Sheet is available on request for professional users.  
A Product Data Sheet is available on request.

**Ecolab Europe GmbH**  
Richtistrasse 7  
CH-8304 Wallisellen  
www.ecolab.com