




HorolithStar



Top performance acidic cleaner for milking equipment and milk tanks.

- ▲ Strong performance in descaling
- ▲ Perfect balance of three acids
- ▲ Best added value solution

APPLICATION

Pre-clean with water. Clean with **HorolithStar** at 0.5% (minimum 10 min. at 40 - 50°C). Rinse with water of drinking quality. To be used in daily rotation with **AseptoStar**.



FEATURES AND BENEFITS

HorolithStar is a well-balanced mixture of three acid components which assures excellent scale removing and cleaning properties. It is suitable for all kind of standard milking equipment, including cooling tanks.

The **Star** products are the top level program to ensure the best hygiene results.

SPECIFICATIONS

HorolithStar	
Physical form	Colourless liquid
Active substance	Inorganic acids
pH at 20°C:	1.0 - 1.5
Density at 20°C (g/cm ³):	1.18 - 1.22
Storage temperature:	-20°C to +40°C

LOGISTICS

Packaging	12 kg	24 kg	72 kg	245 kg
Product code	2332120	2332130	2334760	2332150
EAN code	4028162332124	4028162332131	4028162334760	4028162332155

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtistrasse 7
CH-8304 Wallisellen
www.ecolab.com